



Mouse Anti-human CD4/FITC antibody (BF00201)

BF00201

BF00201

Product Name:	human CD4/FITC
Chinese Name:	FITC标记小鼠抗人CD4单克隆抗体
Organism Species:	Mouse
Clonality:	Monoclonal
克隆号:	?HIT4a
React Species:	Human,Human,
Applications:	Flow Cyt=5ul/Test,Flow-Cyt=5ul/Test not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	56kDa
Form:	Liquid
immunogen:	:
Lsotype:	IgG2b
Storage Buffer:	0.01M PBS, 0.5%BSA, 0.03%Proclin300
Storage:	Store at 2-8°C, protect from light.
PubMed:	PubMed
Product Detail:	<p>In immunology, the CD3 (cluster of differentiation 3) T cell co-receptor helps to activate both the cytotoxic T cell (CD8+ naive T cells) and also T helper cells (CD4+ naive T cells). It consists of a protein complex and is composed of four distinct chains. In mammals, the complex contains a CD3γ chain, a CD3δ chain, and two CD3ϵ chains. These chains associate with the T-cell receptor (TCR) and the ζ-chain (zeta-chain) to generate an activation signal in T lymphocytes. The TCR, ζ-chain, and CD3 molecules together constitute the TCR complex.</p> <p>SWISS: P01730</p> <p>Gene ID: 920</p>

	<p>Important Note:</p>
--	-------------------------------

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com