



Mouse Anti-human CD8/Cy5 antibody (BF00313)

BF00313

BF00313

Product Name:	human CD8/Cy5
Chinese Name:	Cy5标记小鼠抗人CD8单克隆抗体
Organism Species:	Mouse
Clonality:	Monoclonal
克隆号:	?HIT8a
React Species:	
Applications:	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	32-34kDa
Form:	Liquid
immunogen:	:
Lsotype:	IgG1
Storage Buffer:	0.01M PBS, 0.5%BSA, 0.03%Proclin300
Storage:	Store at 2-8°C, protect from light.
PubMed:	PubMed
Product Detail:	<p>In immunology, the CD3 (cluster of differentiation 3) T cell co-receptor helps to activate both the cytotoxic T cell (CD8+ naive T cells) and also T helper cells (CD4+ naive T cells). It consists of a protein complex and is composed of four distinct chains. In mammals, the complex contains a CD3γ chain, a CD3δ chain, and two CD3ϵ chains. These chains associate with the T-cell receptor (TCR) and the ζ-chain (zeta-chain) to generate an activation signal in T lymphocytes. The TCR, ζ-chain, and CD3 molecules together constitute the TCR complex.</p> <p>SWISS: P01732</p> <p>Gene ID: 925</p>

	<p>Important Note:</p>
--	-------------------------------

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com