



Mouse Anti-human CD22/PE-Cy7 antibody (BF00819)

BF00819

BF00819

Product Name:	human CD22/PE-Cy7
Chinese Name:	PE-Cy7标记小鼠抗人CD22单克隆抗体
Organism Species:	Mouse
Clonality:	Monoclonal
克隆号:	?HIB22
React Species:	
Applications:	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	140kDa
Form:	Liquid
immunogen:	:
Lsotype:	?IgG1
Storage Buffer:	0.01M PBS, 0.5%BSA, 0.03%Proclin300
Storage:	Store at 2-8°C, protect from light.
PubMed:	PubMed
Product Detail:	<p>In immunology, the CD3 (cluster of differentiation 3) T cell co-receptor helps to activate both the cytotoxic T cell (CD8+ naive T cells) and also T helper cells (CD4+ naive T cells). It consists of a protein complex and is composed of four distinct chains. In mammals, the complex contains a CD3γ chain, a CD3δ chain, and two CD3ϵ chains. These chains associate with the T-cell receptor (TCR) and the ζ-chain (zeta-chain) to generate an activation signal in T lymphocytes. The TCR, ζ-chain, and CD3 molecules together constitute the TCR complex.</p> <p>SWISS: P20273</p>

Gene ID:
933

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com