



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

[KD-Validated] Anti-NDRG1 Rabbit Monoclonal Antibody

Cat No.: KD-10962

Aliases:

NDRG1; Differentiation-related gene 1 protein; DRG-1; RTP; Rit42; NDRG1; CAP43; DRG1; RTP; targ1; TDD5; tdds

Background:

UniProt Entry: [Q92597](#);NCBI Gene Entry: [10397](#)

Application Information

Molecular Weight: Predicted, 43 kDa; observed, 43 kDa Clonality: Rabbit monoclonal antibody Clone ID: 24GB8605 Species Reactivity: Human, mouse Applications Tested: Western blotting (WB)

Immunogen

A synthesized peptide derived from human NDRG1

Isotype

Rabbit IgG

Storage Buffer

Supplied in PBS (pH 7.4) containing 50% glycerol, and 0.02% sodium azide.

Storage


Store at -20 °C for one year.

Recommended Dilutions

Western Blotting (WB): 1:1,000-1:5,000

Protocols

For general and specific antibody protocols please visit our website. Read all instructions before using this product.



Western blotting analysis using anti-NDRG1 antibody 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQTM ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQTM ECL Substrate Kit 1:5,000 and HRP-conjugated goat anti-rabbit secondary antibody 1:20,000 respectively. Image was developed using NaQTM ECL Substrate Kit .