



Goat Anti-Rabbit IgG/Alexa Fluor 594 antibody

SL0295G-AF594

Product Name:	Goat Anti-Rabbit IgG/Alexa Fluor 594
Chinese Name:	Alexa Fluor 594标记的山羊抗兔IgG
Organism Species:	Goat
Clonality:	Polyclonal
React Species:	
Applications:	IF=1:100-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form:	Lyophilized or Liquid
Concentration:	2mg/1ml
immunogen:	:
Lsotype:	IgG
Purification:	affinity purified by Protein G
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail:	