



Rabbit anti-chicken IgG whole serum antibody

SL0310Rs

Product Name:	Rabbit anti-chicken IgG whole serum
Chinese Name:	兔抗鸡IgG抗血清
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	chk
Applications:	WB=1:10000-50000ELISA=1:10000-100000IHC-P=1:10000-50000IHC-F=1:10000-50000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	180kDa
Form:	lyophilized powder
immunogen:	Chicken IgG purified purified from serum:
Lsotype:	whole antiserum
Purification:	Unpurified
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Product Detail:	