



Rabbit Anti-Hepatitis B Virus Surface Antigen antibody

SL0331R

Product Name:	Hepatitis B Virus Surface Antigen
Chinese Name:	肝炎病毒S蛋白抗体
Alias:	HBSAG_HBVE1; Large envelope protein; L glycoprotein; L-HBsAg; LHB; Large S protein; Large surface protein; Major surface antigen.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Hepatitis B virus (HBVE1)
Applications:	IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	44kDa
Cellular localization:	cytoplasmicThe cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from HBVE1 Hepatitis B Virus Surface Antigen:
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	