



Rabbit Anti-T3 antibody

SL0339R

Product Name:	T3
Chinese Name:	三碘甲状腺素T3抗体
Alias:	3,3',5-Triiodo-L-Thyronine (T3); TRIIODOTHYROXINE.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	T3
Applications:	ELISA=1:500-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	0.65097kDa
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	BSA conjugated T3:T3
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	T3 (Triiodothyroxine) although only a minority of the total thyroxine, but the T3 is very important to check thyroid function. In the biochemical function of T4 than strong, and in the case of hyperthyroidism, T3 concentration in the blood are often marked increase (or even up to 5 times the normal). SWISS: N/A

CAS:
6893-02-3

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com