



Rabbit Anti-Nonylphenol antibody

SL0353R

Product Name:	Nonylphenol
Chinese Name:	壬基酚抗体
Alias:	NP.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Nonylphenol
Applications:	ELISA=1:500-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH-Nonylphenol:
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	药物 化合物抗体 壬基酚是环境荷爾蒙, 現多用在清洁剂中, 做为非离子界面活性剂, 具有类似人类雌激素的效果, 一旦进入人体, 会影响干扰内分泌系統, 危害健康, 甚至造成雄性雌性化。