



## Rabbit Anti-Collagen II antibody

SL0709R

<b>Product Name:</b>	Collagen II
<b>Chinese Name:</b>	II 型胶原 $\alpha$ 1 蛋白/软骨钙素抗体
<b>Alias:</b>	Collagen II alpha 1; COL2A1; COL2A1 protein; collagen, type II, alpha 1; collagen alpha-1(II); type II collagen; alpha-1 type II collagen; alpha1 type II collagen; Col2a1; AOM; Cartilage collagen; Chondrocalcin; COL11A3; Collagen alpha 1(II) chain precursor; Collagen II alpha 1 polypeptide; Collagen type II alpha 1 (primary osteoarthritis spondyloepiphyseal dysplasia congenital); MGC131516; SEDC; Collagen alpha-1(II) chain; Alpha-1 type II collagen; CO2A1 HUMAN.
<b>文献引用</b> <b>PubMed</b> :	<p><b>Specific References(10)</b>SL0709R has been referenced in 10 publications.</p> <p><b>[IF=7.60]</b>Wang, Yang, et al. "Evaluation of three-dimensional scaffolds prepared from poly (3-hydroxybutyrate-co-3-hydroxyhexanoate) for growth of allogeneic chondrocytes for cartilage repair in rabbits." <i>Biomaterials</i> 29.19 (2008): 2858-2868.<b>IHC-P;Rabbit.</b> <a href="#">PubMed:18406457</a></p> <p><b>[IF=3.82]</b>Luo, Yang, et al. "Alendronate Retards the Progression of Lumbar Intervertebral Disc Degeneration in Ovariectomized Rats." <i>Bone</i> (2013).<b>IHC-P;Rat.</b> <a href="#">PubMed:23500174</a></p> <p><b>[IF=5.26]</b>Li, Guo, et al. "Poly (3-hydroxybutyrate-co-4-hydroxybutyrate) Based Electrospun 3D Scaffolds for Delivery of Autogeneic Chondrocytes and Adipose-Derived Stem Cells: Evaluation of Cartilage Defects in Rabbit." <i>Journal of Biomedical Nanotechnology</i> 11.1 (2015): 105-116.<b>other;Rabbit.</b> <a href="#">PubMed:26301304</a></p> <p><b>[IF=3.71]</b>Fu, Na, et al. "P34HB film promotes cell adhesion, in vitro proliferation, and in</p>

	<p>vivo cartilage repair." RSC Advances (2015).<b>IHC-P;Rat.</b></p> <p><a href="#">PubMed:not posted yet</a></p> <p><b>[IF=4.25]</b>Chen, Ying-Chen, et al. "Development and characterization of acellular extracellular matrix scaffolds from porcine menisci for use in cartilage tissue engineering." Tissue Engineering ja (2015).<b>other;</b></p> <p><a href="#">PubMed:25919905</a></p> <p><b>[IF=1.27]</b>Liu, Yong, et al. "Lentivirus-mediated TGF-β3, CTGF and TIMP1 gene transduction as a gene therapy for intervertebral disc degeneration in an in vivo rabbit model." Experimental and Therapeutic Medicine.<b>WB;Rabbit.</b></p> <p><a href="#">PubMed:27073456</a></p> <p><b>[IF=0.00]</b>Sheu, Ming-thau, and Ying-chen Chen. "METHOD FOR PRODUCTION OF DECELLULARIZED BIOLOGICAL MATERIAL AND THE DECELLULARIZED BIOLOGICAL MATERIAL PREPARED THEREFROM." U.S. Patent No. 20,150,344,842. 3 Dec. 2015.<b>IHC-P;Pig.</b></p> <p><a href="#">PubMed:1111</a></p> <p><b>[IF=4.45]</b>Wang, Chen-Chie, et al. "Expandable Scaffold Improves Integration of Tissue Engineered Cartilage: An In Vivo Study in A Rabbit Model." Tissue Engineering ja (2016).<b>IHC-P;Rabbit.</b></p> <p><a href="#">PubMed:27193498</a></p> <p><b>[IF=1.50]</b>Zheng, Li, et al. "Comparative study of collagen hydrogels modified in two ways using the model of ectopic cartilage construction with diffusion-chamber in immunocompetent host." Journal of applied biomaterials &amp; functional materials 12.1 (2014).<b>IHC-P;Rat.</b></p> <p><a href="#">PubMed:22865571</a></p> <p><b>[IF=3.08]</b>Fu, Na, et al. "PCL-PEG-PCL film promotes cartilage regeneration in vivo."Cell Proliferation (2016).<b>IHC-P;Rat.</b></p> <p><a href="#">PubMed:27647680</a></p>
<b>Organism Species:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>React Species:</b>	Human,Mouse,Rat,Chicken,Dog,Pig,Cow,Rabbit,Guinea Pig,
<b>Applications:</b>	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

<b>Molecular weight:</b>	117kDa
<b>Cellular localization:</b>	Extracellular matrixSecretory protein
<b>Form:</b>	Lyophilized or Liquid
<b>Concentration:</b>	1mg/ml
<b>immunogen:</b>	KLH conjugated synthetic peptide derived from human Collagen II:1201-1300/1487
<b>Lsotype:</b>	IgG
<b>Purification:</b>	affinity purified by Protein A
<b>Storage Buffer:</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Storage:</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>PubMed:</b>	<a href="#">PubMed</a>
<b>Product Detail:</b>	<p>This gene encodes the alpha-1 chain of type II collagen, a fibrillar collagen found in cartilage and the vitreous humor of the eye. Mutations in this gene are associated with achondrogenesis, chondrodysplasia, early onset familial osteoarthritis, SED congenita, Langer-Saldino achondrogenesis, Kniest dysplasia, Stickler syndrome type I, and spondyloepimetaphyseal dysplasia Strudwick type. In addition, defects in processing chondrocalcin, a calcium binding protein that is the C-propeptide of this collagen molecule, are also associated with chondrodysplasia. There are two transcripts identified for this gene. [provided by RefSeq, Jul 2008]</p> <p><b>Function:</b> Type II collagen is specific for cartilaginous tissues. It is essential for the normal embryonic development of the skeleton, for linear growth and for the ability of cartilage to resist compressive forces.</p> <p><b>Subunit:</b> Homotrimers of alpha 1(II) chains.</p> <p><b>Subcellular Location:</b> Secreted, extracellular space, extracellular matrix.</p> <p><b>Tissue Specificity:</b> Isoform 2 is highly expressed in juvenile chondrocyte and low in fetal chondrocyte.</p> <p><b>Post-translational modifications:</b> Proline residues at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains. Proline residues at the second position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some of the chains. The N-telopeptide is covalently linked to the helical COL2 region of alpha 1(IX), alpha 2(IX) and alpha 3(IX) chain. The C-telopeptide is covalently linked to an another site in the helical region of alpha 3(IX) COL2.</p> <p><b>DISEASE:</b> Defects in COL2A1 are the cause of spondyloepiphyseal dysplasia congenital type</p>

(SEDC) [MIM:183900]. This disorder is characterized by disproportionate short stature and pleiotropic involvement of the skeletal and ocular systems.

Defects in COL2A1 are the cause of spondyloepimetaphyseal dysplasia Strudwick type (SEMD-STR) [MIM:184250]. A bone disease characterized by disproportionate short stature from birth, with a very short trunk and shortened limbs, and skeletal abnormalities including lordosis, scoliosis, flattened vertebrae, pectus carinatum, coxa vara, clubfoot, and abnormal epiphyses or metaphyses. A distinctive radiographic feature is irregular sclerotic changes, described as dappled in the metaphyses of the long bones.

Defects in COL2A1 are the cause of achondrogenesis type 2 (ACG2) [MIM:200610]; also known as achondrogenesis-hypochondrogenesis type II. ACG2 is a disease characterized by the absence of ossification in the vertebral column, sacrum and pubic bones.

Defects in COL2A1 are the cause of Legg-Calve-Perthes disease (LCPD) [MIM:150600]; also known as Legg-Perthes disease or Perthes disease. LCPD is characterized by loss of circulation to the femoral head, resulting in avascular necrosis in a growing child. Clinical pictures of the disease vary, depending on the phase of disease progression through ischemia, revascularization, fracture and collapse, and repair and remodeling of the bone.

Defects in COL2A1 are the cause of Kniest dysplasia (KD) [MIM:156550]; also known as Kniest syndrome or metatropic dwarfism type II. KD is a moderately severe chondrodysplasia phenotype that results from mutations in the COL2A1 gene. Characteristics of the disorder include a short trunk and extremities, mid-face hypoplasia, cleft palate, myopia, retinal detachment, and hearing loss.

Defects in COL2A1 are a cause of primary avascular necrosis of femoral head (ANFH) [MIM:608805]; also known as ischemic necrosis of the femoral head or osteonecrosis of the femoral head. ANFH causes disability that often requires surgical intervention. Most cases are sporadic, but families in which there is an autosomal dominant inheritance of the disease have been identified. It has been estimated that 300,000 to 600,000 people in the United States have ANFH. Approximately 15,000 new cases of this common and disabling disorder are reported annually. The age at the onset is earlier than that for osteoarthritis. The diagnosis is typically made when patients are between the ages of 30 and 60 years. The clinical manifestations, such as pain on exertion, a limping gait, and a discrepancy in leg length, cause considerable disability. Moreover, nearly 10 percent of the 500,000 total-hip arthroplasties performed each year in the United States involve patients with ANFH. As a result, this disease creates a substantial socioeconomic cost as well as a burden for patients and their families.

Defects in COL2A1 are the cause of osteoarthritis with mild chondrodysplasia (OACD) [MIM:604864]. Osteoarthritis is a common disease that produces joint pain and stiffness together with radiologic evidence of progressive degeneration of joint cartilage. Some forms of osteoarthritis are secondary to events such as trauma, infections, metabolic disorders, or congenital or heritable conditions that deform the epiphyses or related structures. In most patients, however, there is no readily identifiable cause of osteoarthritis. Inheritance in a Mendelian dominant manner has been demonstrated in some families with primary generalized osteoarthritis. Reports demonstrate coinheritance of primary generalized osteoarthritis with specific alleles of the gene COL2A1, the precursor of the major protein of cartilage.

Defects in COL2A1 are the cause of platyspondylic lethal skeletal dysplasia Torrance type (PLSD-T) [MIM:151210]. Platyspondylic lethal skeletal dysplasias (PLSDs) are a

heterogeneous group of chondrodysplasias characterized by severe platyspondyly and limb shortening. PLSD-T is characterized by varying platyspondyly, short ribs with anterior cupping, hypoplasia of the lower ilia with broad ischial and pubic bones, and shortening of the tubular bones with splayed and cupped metaphyses. Histology of the growth plate typically shows focal hypercellularity with slightly enlarged chondrocytes in the resting cartilage and relatively well-preserved columnar formation and ossification at the chondro-osseous junction. PLSD-T is generally a perinatally lethal disease, but a few long-term survivors have been reported.

Defects in COL2A1 are the cause of multiple epiphyseal dysplasia with myopia and conductive deafness (EDMMD) [MIM:132450]. Multiple epiphyseal dysplasia is a generalized skeletal dysplasia associated with significant morbidity. Joint pain, joint deformity, waddling gait, and short stature are the main clinical signs and symptoms. EDMMD is an autosomal dominant disorder characterized by epiphyseal dysplasia associated with progressive myopia, retinal thinning, crenated cataracts, conductive deafness.

Defects in COL2A1 are the cause of spondyloperipheral dysplasia (SPD) [MIM:271700]. SPD patients manifest short stature, midface hypoplasia, sensorineural hearing loss, spondyloepiphyseal dysplasia, platyspondyly and brachydactyly.

**Similarity:**

Belongs to the fibrillar collagen family.

Contains 1 fibrillar collagen NC1 domain.

Contains 1 VWFC domain.

**SWISS:**

P02458

**Gene ID:**

1280

**Database links:**

[Entrez Gene: 1280](#) Human

[Entrez Gene: 12824](#) Mouse

[Entrez Gene: 25412](#) Rat

[Omim: 120140](#) Human

[SwissProt: P02458](#) Human

[SwissProt: P28481](#) Mouse

[SwissProt: P05539](#) Rat

[Unigene: 408182](#) Human

[Unigene: 2423](#) Mouse

[Unigene: 10124](#) Rat

**Important Note:**


This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

II型胶原是软骨基质中的一种结构蛋白；II型胶原不与I、III、IV、V、VI胶原以及其它血清蛋白或非胶原性细胞外相关蛋白起React

Species。II型胶原最先发现于软骨基质中，在眼睛中也有少量存在。组成II型胶原的纤维较I型胶原更纤细。该抗体可以特异性识别II型胶原，与其他类型的胶原无React Species, 抗II型胶原抗体-

主要用于II型胶原分布及变态反应方面及良/恶性Tumour的Extracellular matrix方面的研究。

Picture:



Sample:

293T (human)Cell Lysate at 40 ug

A549 (human)Cell Lysate at 40 ug

MCF-7 (human)Cell Lysate at 40 ug

HepG2 (human)Cell Lysate at 40 ug

Placenta (mouse) Lysate at 40 ug


MG63 (human) Cell Lysate at 40 ug

Primary: Anti- Collagen II (SL0709R) at 1/300 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 117 kD

Observed band size: 120 kD




Tissue/cell: human esophageal carcinoma; 4% Paraformaldehyde-fixed and paraffin-embedded;

Antigen retrieval: citrate buffer ( 0.01M, pH 6.0 ), Boiling bathing for 15min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30min; Blocking buffer (normal goat serum,C-0005) at 37°C for 20 min;

Incubation: Anti-Collagen II Polyclonal Antibody, Unconjugated(SL0709R) 1:200, overnight at 4°C, followed by conjugation to the secondary antibody(SP-0023) and

## DAB(C-0010) staining



Tissue/cell: rabbit meniscus tissue; 4% Paraformaldehyde-fixed and paraffin-embedded;

Antigen retrieval: citrate buffer ( 0.01M, pH 6.0 ), Boiling bathing for 15min;

Blocking buffer (normal goat serum, C-0005) at 37°C for 20 min;

Incubation: Anti-Collagen II Polyclonal Antibody, Unconjugated(SL0709R) 1:200, overnight at 4°C; The secondary antibody was Goat Anti-Rabbit IgG, Cy3 conjugated(SL0709R) used at 1:200 dilution for 40 minutes at 37°C.

DAPI(5ug/ml, blue, C-0033) was used to stain the cell nuclei