



Rabbit Anti-Os05g0477200 protein antibody

SL0807R

| | |
|--------------------------|---|
| Product Name: | Os05g0477200 protein |
| Chinese Name: | 水稻白叶枯病菌Os05g0477200蛋白抗体 |
| Alias: | Os05g0477200 protein; Os05g0477200 [Oryza sativa Japonica Group]. |
| Organism Species: | Rabbit |
| Clonality: | Polyclonal |
| React Species: | Oryza sativa Japonica Group |
| Applications: | WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Molecular weight: | 52kDa |
| Form: | Lyophilized or Liquid |
| Concentration: | 1mg/ml |
| immunogen: | KLH conjugated synthetic peptide derived from the middle of Oryza sativa Os05g0477200 protein:221-320/472 |
| Lsotype: | IgG |
| Purification: | affinity purified by Protein A |
| Storage Buffer: | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage: | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| PubMed: | PubMed |
| Product Detail: | |