



## Rabbit Anti-HHV8 ORF50 antibody

SL0860R

<b>Product Name:</b>	HHV8 ORF50
<b>Chinese Name:</b>	人类疱疹病毒8抗体
<b>Alias:</b>	HHV8 ORF50; HHV 8; Human Herpes Virus 8; Human herpesvirus 8; Kaposi's sarcoma associated herpes virus; Kaposi's sarcoma associated herpesvirus; KSHV.
<b>Organism Species:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>React Species:</b>	herpesvirus8typeP
<b>Applications:</b>	ELISA=1:500-1000IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight:</b>	74kDa
<b>Form:</b>	Lyophilized or Liquid
<b>Concentration:</b>	1mg/ml
<b>immunogen:</b>	KLH conjugated synthetic peptide derived from HHV8/ORF50:601-691/691
<b>Lsotype:</b>	IgG
<b>Purification:</b>	affinity purified by Protein A
<b>Storage Buffer:</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Storage:</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>PubMed:</b>	<a href="#">PubMed</a>
<b>Product Detail:</b>	HHV8 has been found to be associated with three different diseases observed in AIDS patients; kaposi's sarcoma, primary effusion lymphoma (which is a rare type of non-Hodgkin lymphoma affecting the body cavities) and multicentric Castleman's disease. To date there is much evidence to support a direct role for HHV8 in kaposi's sarcoma.  <b>SWISS:</b> N/A

**Gene ID:**

N/A

**Database links:**

NCBI Reference Sequence: YP\_001129401

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

学者研究认为:包括HHV-

8在内的一类疱疹病毒THC(Herpesviruses)都存在ORF50基因的诱导, ORF50能够靶向这类病毒。

www.sunlongbiotech.com