



Rabbit Anti-PRRSV antibody

SL10043R

Product Name:	PRRSV
Chinese Name:	猪蓝耳病病毒抗体
Alias:	Porcine reproductive and respiratory syndrome, PRRS.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	PRRSV
Applications:	ELISA=1:500-1000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	PRRSV vaccine:
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	猪繁殖与呼吸综合征(Porcine reproductive and respiratory syndrome, PRRS)也称猪蓝耳病, 是由动脉炎病毒科动脉炎病毒属的猪繁殖与呼吸综合征病毒,(Porcine reproductive and respiratory syndrome virus, PRRSV)引起的一种严重的病毒性传染病, 该病以引起母猪的繁殖障碍、仔猪和育成猪呼吸道症状及高死亡率为主要特征。自1987年首次发现以来, 该病广泛流行于世界各主要养猪国家和地区, 已成为威胁养猪业的主要疫病之一。2006年以来, 由该病毒变异株引起的高致病性猪蓝耳病(Highly pathogenic porcine reproductive and respiratory syndrome, HPPRRS)在我国广泛流行给养猪业带来了巨大的损失。

