



Rabbit Anti-Interferon alpha 8 antibody

SL10573R

| | |
|-------------------------------|---|
| Product Name: | Interferon alpha 8 |
| Chinese Name: | Interferon α 8抗体 |
| Alias: | IFN Alpha 8; IFNA8; IFN-alpha-8; IFN-alpha 8; IFNA8; IFN-alphaB; Interferon alpha B; Interferon alpha B2; interferon alpha type 201; interferon alpha-B'; interferon, alpha 8; LeIF B; OTTHUMP00000021146; IFNA8_HUMAN. |
| Organism Species: | Rabbit |
| Clonality: | Polyclonal |
| React Species: | Human, |
| Applications: | ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Molecular weight: | 18kDa |
| Cellular localization: | Secretory protein |
| Form: | Lyophilized or Liquid |
| Concentration: | 1mg/ml |
| immunogen: | KLH conjugated synthetic peptide derived from human Interferon alpha 8: |
| Lsotype: | IgG |
| Purification: | affinity purified by Protein A |
| Storage Buffer: | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage: | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| PubMed: | PubMed |
| Product Detail: | |