



Rabbit Anti-PRRSV N protein antibody

SL10766R

Product Name:	PRRSV N protein
Chinese Name:	猪蓝耳病病毒N蛋白抗体
Alias:	N [Porcine reproductive and respiratory syndrome virus]; nucleoprotein [Porcine reproductive and respiratory syndrome virus]; PRRSV nucleoprotein.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	PRRSV
Applications:	ELISA=1:5000-10000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	13kDa
Cellular localization:	The nucleus
Form:	Liquid
Concentration:	1mg/ml
immunogen:	Recombinant PRRSV N protein:
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	PRRSV is a small, enveloped RNA virus. It contains a single-stranded, positive-sense, RNA genome with a size of approximately 15 kilobases. The genome contains nine open reading frames. PRRSV is a member of the genus Arterivirus, family Arteriviridae, order Nidovirales. Subclinical infections are common, with clinical signs occurring sporadically in a herd. Clinical signs include reproductive failure in sows such as abortions and giving birth to stillborn or mummified fetuses, and cyanosis of the ear and vulva. In neonatal pigs, the disease causes respiratory distress, with

increased susceptibility to respiratory infections such as Glasser's disease.

SWISS:

N/A

Gene ID:

N/A

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com