



Rabbit Anti-AADACL2 antibody

SL11385R

Product Name:	AADACL2
Chinese Name:	芳香乙酰胺脱乙酰基酶样蛋白2抗体
Alias:	AADACL2; ADCL2_HUMAN; Arylacetamide deacetylase like 2 [Precursor]; Arylacetamide deacetylase-like 2; MGC72001.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Cow,Rabbit,Sheep,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	44kDa
Cellular localization:	Secretory protein
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human AADACL2:71-115/401
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	No data available. Subcellular Location: Secreted. Similarity:

Belongs to the 'GDXG' lipolytic enzyme family.

SWISS:

Q6P093

Gene ID:

344752

Database links:

[Entrez Gene: 344752](#)Human

[GenBank: NM_207365](#)Human

[SwissProt: Q6P093](#)Human

[Unigene: 144710](#)Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com