



Rabbit Anti-Nogo B receptor antibody

SL11468R

Product Name:	Nogo B receptor
Chinese Name:	轴索过度生长抑制因子B受体抗体
Alias:	C6orf68; NGBR; nuclear undecaprenyl pyrophosphate synthase 1 homolog; NUS1; NGBR_HUMAN.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Dog,Pig,Cow,Horse,Sheep,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	30kDa
Cellular localization:	cytoplasmicThe cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human NGBR/Nogo B receptor:167-210/293
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	Nogo is an oligodendrocyte-specific member of the Reticulon family and is a component of CNS white matter that inhibits axon outgrowth, induces collapse of growth cones of chick dorsal root ganglion cells, and inhibits the spreading of 3T3 fibroblasts. Nogo is expressed by oligodendrocytes but not by Schwann cells and associates primarily with the endoplasmic reticulum. Nogo exists in three different

splice forms, Nogo-A, -B and -C. NgBR (Nogo-B receptor), also known as nuclear undecaprenyl pyrophosphate synthase 1 homolog, is a 293 amino acid single-pass type I membrane protein that acts as a specific receptor for the amino-terminus of Nogo-B. Through this interaction, NgBR is involved in the regulation of vascular remodeling and angiogenesis. NgBR also enhances Niemann-Pick type C2 protein (NPC2) stabilization. Knockdown of NgBR mRNA leads to decreased NPC2 levels, which results in the hallmarks of NPC2 mutation: increased intracellular cholesterol accumulation and a loss of sterol sensing.

Function:

Acts as a specific receptor for the N-terminus of Nogo-B, a neural and cardiovascular regulator. Able to regulate vascular remodeling and angiogenesis. Its similarity with UPP synthetase proteins suggests that it may act as a scaffold for the binding of isoprenyl lipids and/or prenylated proteins.

Subunit:

Interacts with DHDDS, promoting its isoprenyltransferase activity. Interacts with NPC2.

Subcellular Location:

Endoplasmic reticulum membrane; Single-pass type I membrane protein.

Note=Colocalizes with Nogo-B during VEGF and wound healing angiogenesis.

Similarity:

Belongs to the UPP synthase family.

SWISS:

Q96E22

Gene ID:

116150

Database links:

[Entrez Gene: 116150](#) Human

[Entrez Gene: 52014](#) Mouse

[Entrez Gene: 294400](#) Rat

[Omim: 610463](#) Human

[SwissProt: Q96E22](#) Human

[SwissProt: Q99LJ8](#) Mouse

[Unigene: 289008](#) Human

[Unigene: 525826](#) Human

[Unigene: 199964](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com