



Rabbit Anti-phospho-ZNF598 (Tyr306) antibody

SL12222R

Product Name:	phospho-ZNF598 (Tyr306)
Chinese Name:	磷酸化Zinc finger protein598抗体
Alias:	ZNF598 (phospho Y306); p-ZNF598 (phospho Y306); p-ZNF598 (phospho Y306); Zinc finger protein 598; ZN598 HUMAN; Znf598.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Horse,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	99kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthesised phosphopeptide derived from human ZNF598 around the phosphorylation site of Tyr306:ED(P-Y)EE
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This protein and Grb10-interacting GYF protein 2 have been identified as a components of the mammalian 4EHP (m4EHP) complex. The complex is thought to function as a

translation repressor in embryonic development. [provided by RefSeq, Oct 2012].

Similarity:

Similarity : Contains 1 C2H2-type zinc finger.

Contains 1 RING-type zinc finger.

SWISS:

Q86UK7

Gene ID:

90850

Database links:

[Entrez Gene: 90850](#) Human

[Entrez Gene: 213753](#) Mouse

[SwissProt: Q86UK7](#) Human

[SwissProt: Q80YR4](#) Mouse

[Unigene: 343828](#) Human

[Unigene: 219581](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.