



Rabbit Anti-phospho-AQP2 (Ser256) antibody

SL12507R

Product Name:	phospho-AQP2 (Ser256)
Chinese Name:	磷酸化水Channel protein2抗体
Alias:	Aquaporin 2 (phospho S256); Aquaporin 2 (phospho Ser256); p-Aquaporin 2 (phospho S256); p-Aquaporin 2 (S256); ADH water channel; AQP 2; AQP CD; AQP2; AQPCD; Aquaporin 2 collecting duct; Aquaporin CD; Aquaporin2; Aquaporine 2; Collecting duct water channel protein; MGC34501; Water channel protein for renal collecting duct; WCH CD; WCHCD; AQP2_HUMAN.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit, Monkey,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	29kDa
Cellular localization:	cytoplasmicThe cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthesised phosphopeptide derived from human AQP2 around the phosphorylation site of Ser256:RQ(p-S)VE
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	Aquaporins (AQPs) are a large family of integral membrane water transport channel

proteins that facilitate the transport of water through the cell membrane. This function is conserved in animals, plants and bacteria. Many isoforms of aquaporin have been identified in mammals, designated AQP0 through AQP10. Aquaporins are widely distributed and it is not uncommon for more than one type of AQP to be present in the same cell. Although most aquaporins are only permeable to water, AQP3, AQP7, AQP9 and one of the two AQP10 transcripts are also permeable to urea and glycerol. AQP2 is the only water channel that is activated by vasopressin to enhance water reabsorption in the kidney collecting duct. Aquaporins are involved in renal water absorption, generation of pulmonary secretions, lacrimation, and the secretion and reabsorption of cerebrospinal fluid and aqueous humor.

Function:

Forms a water-specific channel that provides the plasma membranes of renal collecting duct with high permeability to water, thereby permitting water to move in the direction of an osmotic gradient.

Subcellular Location:

Apical cell membrane. Cytoplasmic vesicle membrane. Shuttles from vesicles to the apical membrane.

Tissue Specificity:

Expressed in renal collecting tubules.

Post-translational modifications:

Ser-256 phosphorylation is necessary and sufficient for expression at the apical membrane. Endocytosis is not phosphorylation-dependent.

DISEASE:

Defects in AQP2 are the cause of diabetes insipidus nephrogenic autosomal (ANDI) [MIM:125800]; also known as diabetes insipidus nephrogenic type 2. ANDI is caused by the inability of the renal collecting ducts to absorb water in response to arginine vasopressin. It is characterized by excessive water drinking (polydypsia), excessive urine excretion (polyuria), persistent hypotonic urine, and hypokalemia. Inheritance can be autosomal dominant or recessive.

Similarity:

Belongs to the MIP/aquaporin (TC 1.A.8) family.

SWISS:

P41181

Gene ID:

359


Database links:

[Entrez Gene: 359](#) Human
[Entrez Gene: 11827](#) Mouse
[Entrez Gene: 25386](#) Rat
[Omim: 107777](#) Human
[SwissProt: P41181](#) Human
[SwissProt: P56402](#) Mouse
[SwissProt: P34080](#) Rat
[Unigene: 130730](#) Human
[Unigene: 20206](#) Mouse
[Unigene: 90076](#) Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Picture:



Sample:

MDCK(Dog) Cell Lysate at 30 ug

Primary: Anti-phospho-AQP2 (Ser256) (SL12507R) at 1/500 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 29 kD

Observed band size: 29 kD

www.sunlongbiotech.com