



Rabbit Anti-PDHA2 antibody

SL12596R

Product Name:	PDHA2
Chinese Name:	PDHA2蛋白抗体
Alias:	EC 1.2.4.1; MGC149517; MGC149518; mitochondrial; ODPAT_HUMAN; Pdha2; PDHAL; PDHE1 A type II; PDHE1-A type II; Pyruvate dehydrogenase (lipoamide) alpha 2; Pyruvate dehydrogenase E1 component subunit alpha; Pyruvate dehydrogenase E1 component subunit alpha, testis specific form, mitochondrial; Pyruvate dehydrogenase, E1 alpha polypeptide, testis spec; testis-specific form.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	40kDa
Cellular localization:	cytoplasmic
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human PDHA2:61-160/388
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	