



Rabbit Anti-PDXDC1 antibody

SL12606R

Product Name:	PDXDC1
Chinese Name:	吡哆醛依赖脱羧酶结构域PDXDC1抗体
Alias:	2210010A19Rik; AA415817; KIAA0251; Kiaa0251 hp; LP8165; MGC189517; MGC8182; PDXD1_HUMAN; PDXDC 1; Pdxdc1; Pyridoxal dependent decarboxylase domain containing 1; Pyridoxal-dependent decarboxylase domain-containing protein 1; rCG 49895; RGD1562597.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Dog,Cow,Horse,Rabbit,Sheep,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	87kDa
Cellular localization:	cytoplasmic
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human PDXDC1:651-750/788
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	