



Rabbit Anti-TIE1 antibody

SL1334R

Product Name:	TIE1
Chinese Name:	血管生成素受体1抗体
Alias:	TIE1_HUMAN; JKT 14; JTK14; TIE; Tyrosine Kinase with Immunoglobulin and Epidermal Growth Factor Homology Domains; Tyrosine protein kinase receptor Tie 1; Tyrosine-protein kinase receptor Tie-1; tyrosine kinase receptor 1; tyrosine kinase with immunoglobulin-like and EGF-like domains 1.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	123kDa
Cellular localization:	cytoplasmic
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human Tie1:851-1000/1138
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	TIE1/TIE (tyrosine kinase with Ig and EGF homology domains 1) and TIE2/Tek define a new class of the receptor tyrosine kinase (RTK) subfamily with unique structural characteristics: two immunoglobulin like domains flanking three epidermal growth factor (EGF) like domains followed by three fibronectin type III like repeats in the

extracellular region and a split tyrosine kinase domain in the cytoplasmic region. Human TIE1 cDNA encodes a 1138 amino acid residue precursor protein with a putative signal peptide, an extracellular domain, and a cytoplasmic domain. Human TIE1/Fc, a disulfide linked homodimeric protein, has a calculated molecular mass of approximately 107 kDa. Due to glycosylation, the protein migrates to approximately 160 kDa in SDS PAGE under reducing conditions. TIE1 and TIE2, expressed primarily on endothelial and hematopoietic progenitor cells, play important roles in angiogenesis, vasculogenesis, and hematopoiesis. In developing vascular endothelial cells, TIE1 and TIE2 are specifically expressed. Two ligands that bind TIE have been identified, angiopoietin 1 and angiopoietin 2. Based on gene targeting studies, the in vivo functions of TIE1 are related to endothelial cell differentiation. The receptor tyrosine kinase TIE also plays a role in the survival and integrity of the endothelium.

Function:

Transmembrane tyrosine-protein kinase that may modulate TEK/TIE2 activity and contribute to the regulation of angiogenesis.

Subunit:

Heterodimer with TEK/TIE2.

Subcellular Location:

Cell membrane; Single-pass type I membrane protein.

Tissue Specificity:

Specifically expressed in developing vascular endothelial cells.

Post-translational modifications:

Phosphorylated on tyrosine residues in response to ANGPT1, most likely by TEK/TIE2.

Similarity:

Belongs to the protein kinase superfamily. Tyr protein kinase family. Tie subfamily.

Contains 3 EGF-like domains.

Contains 3 fibronectin type-III domains.

Contains 2 Ig-like C2-type (immunoglobulin-like) domains.

Contains 1 protein kinase domain.

SWISS:

P35590

Gene ID:

7075

Database links:

[Entrez Gene: 280941](#)Cow

[Entrez Gene: 7075](#)Human

[Entrez Gene: 21846](#)Mouse

[Entrez Gene: 89806](#)Rat

[Omim: 600222](#)Human

[SwissProt: Q06805](#)Cow

[SwissProt: P35590](#)Human

[SwissProt: Q06806](#)Mouse

[Unigene: 78824](#)Human

[Unigene: 4345](#)Mouse

[Unigene: 13171](#)Rat


Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

血管生成素-

1受体又称上皮生长因子样域酪氨酸激酶1或内皮酪氨酸激酶受体,Tie1在vascular endothelial cell、腺上皮及间质细胞中均有表达,主要定位于细胞质。

Picture:



Sample:

A549(Human) Cell Lysate at 30 ug


HepG2(Human) Cell Lysate at 30 ug

Primary: Anti-TIE1 (SL1334R) at 1/1000 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 123 kD

Observed band size: 123 kD



Tissue/cell: human brain tissue; 4% Paraformaldehyde-fixed and paraffin-embedded;

Antigen retrieval: citrate buffer (0.01M, pH 6.0), Boiling bathing for 15min; Block endogenous peroxidase by 3% Hydrogen peroxide for 30min; Blocking buffer (normal goat serum,C-0005) at 37°C for 20 min;

Incubation: Anti-TIE1 Polyclonal Antibody, Unconjugated(SL1334R) 1:200, overnight at 4°C, followed by conjugation to the secondary antibody(SP-0023) and DAB(C-0010) staining