



## Rabbit Anti-CSMD1 antibody

SL1579R

<b>Product Name:</b>	CSMD1
<b>Chinese Name:</b>	CSMD1抗体
<b>Alias:</b>	CUB and Sushi multiple domains1; CUB and sushi domain-containing protein 1; CUB and sushi multiple domains protein 1; B930082H09; CSMD1; KIAA1890; Mdcrl; mKIAA1890; CSMD1_HUMAN.
<b>Organism Species:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>React Species:</b>	Human,Mouse,Rat,Chicken,Cow,Horse,
<b>Applications:</b>	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight:</b>	392kDa
<b>Cellular localization:</b>	The cell membrane
<b>Form:</b>	Lyophilized or Liquid
<b>Concentration:</b>	1mg/ml
<b>immunogen:</b>	KLH conjugated synthetic peptide derived from human CSMD1:3510-3565/3566
<b>Lsotype:</b>	IgG
<b>Purification:</b>	affinity purified by Protein A
<b>Storage Buffer:</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Storage:</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>PubMed:</b>	<a href="#">PubMed</a>
<b>Product Detail:</b>	CSMD1基因是一种Transmembrane protein,CSMD1蛋白的下降与多种Tumour的发生密切相关,近年来CSMD1在头颈部Tumour的研究方面较多,成为关注的焦点。