



Rabbit Anti-ZNF77 antibody

SL16482R

Product Name:	ZNF77
Chinese Name:	Zinc finger protein77抗体
Alias:	pT1; Zinc finger protein 77; ZNF77 HUMAN; ZNFpT1.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Horse,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	62kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human ZNF77:451-545/545
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	ZNF77 belongs to the krueppel C2H2 type zinc finger protein family. It contains 12 C2H2 type zinc fingers and 1 KRAB domain. ZNF77 may be involved in transcriptional regulation. Function: May be involved in transcriptional regulation.

Subcellular Location:

Nuclear

Similarity:

Belongs to the krueppel C2H2-type zinc-finger protein family.

Contains 12 C2H2-type zinc fingers.

Contains 1 KRAB domain.

SWISS:

Q15935

Gene ID:

58492

Database links:

[Entrez Gene: 58492](#) Human

[GenBank: BC062561](#) Human

[Omim: 194551](#) Human

[SwissProt: Q15935](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.