



Rabbit Anti-Influenza A PB1-F2 antibody

SL16648R

Product Name:	Influenza A PB1-F2
Chinese Name:	A型流感病毒PB1-F2抗体
Alias:	PB1-F2 protein; PB1F2 I07A0; Protein PB1-F2.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Influenza A
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	6.6kDa
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from Influenza A Influenza A PB1:1-57/57
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 癢 for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20癢. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 癢.
PubMed:	PubMed
Product Detail:	