



## Rabbit Anti-SPAG7 antibody

SL17619R

|                               |                                                                                                                                                                                                                                                                                                                               |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | SPAG7                                                                                                                                                                                                                                                                                                                         |
| <b>Chinese Name:</b>          | 精子相关抗原7抗体                                                                                                                                                                                                                                                                                                                     |
| <b>Alias:</b>                 | ACRP; FSA 1; FSA1; Spag7; SPAG7 HUMAN; Sperm-associated antigen 7.                                                                                                                                                                                                                                                            |
| <b>Organism Species:</b>      | Rabbit                                                                                                                                                                                                                                                                                                                        |
| <b>Clonality:</b>             | Polyclonal                                                                                                                                                                                                                                                                                                                    |
| <b>React Species:</b>         | Human,Mouse,Rat,Dog,Pig,Cow,Horse,Rabbit,Sheep,                                                                                                                                                                                                                                                                               |
| <b>Applications:</b>          | ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair)<br>not yet tested in other applications.<br>optimal dilutions/concentrations should be determined by the end user.                                                                                            |
| <b>Molecular weight:</b>      | 26kDa                                                                                                                                                                                                                                                                                                                         |
| <b>Cellular localization:</b> | The nucleus                                                                                                                                                                                                                                                                                                                   |
| <b>Form:</b>                  | Lyophilized or Liquid                                                                                                                                                                                                                                                                                                         |
| <b>Concentration:</b>         | 1mg/ml                                                                                                                                                                                                                                                                                                                        |
| <b>immunogen:</b>             | KLH conjugated synthetic peptide derived from human SPAG7:151-227/227                                                                                                                                                                                                                                                         |
| <b>Lsotype:</b>               | IgG                                                                                                                                                                                                                                                                                                                           |
| <b>Purification:</b>          | affinity purified by Protein A                                                                                                                                                                                                                                                                                                |
| <b>Storage Buffer:</b>        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                                                                                                                                                                                                                                                              |
| <b>Storage:</b>               | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| <b>PubMed:</b>                | <a href="#">PubMed</a>                                                                                                                                                                                                                                                                                                        |
| <b>Product Detail:</b>        |                                                                                                                                                                                                                                                                                                                               |