



Rabbit Anti-SPIN3 antibody

SL17667R

| | |
|-------------------------------|---|
| Product Name: | SPIN3 |
| Chinese Name: | 纺锤体样蛋白3抗体 |
| Alias: | SPIN-3; SPIN3; SPIN3 HUMAN; Spindlin-3; Spindlin-like protein 3. |
| Organism Species: | Rabbit |
| Clonality: | Polyclonal |
| React Species: | Human,Pig, |
| Applications: | ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Molecular weight: | 29kDa |
| Cellular localization: | The nucleus |
| Form: | Lyophilized or Liquid |
| Concentration: | 1mg/ml |
| immunogen: | KLH conjugated synthetic peptide derived from human SPIN3:1-100/258 |
| Lsotype: | IgG |
| Purification: | affinity purified by Protein A |
| Storage Buffer: | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage: | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| PubMed: | PubMed |
| Product Detail: | |