



Rabbit Anti-OR5W2 antibody

SL17955R

Product Name:	OR5W2
Chinese Name:	嗅觉受体5W2抗体
Alias:	Olfactory receptor 5W2; Olfactory receptor 5W3; Olfactory receptor OR11-155; Olfactory receptor, family 5, subfamily W, member 2; Olfactory receptor, family 5, subfamily W, member 3 pseudogene; OR5W2; OR5W2_HUMAN; OR5W2P; OR5W3P.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	35kDa
Cellular localization:	The cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human OR5W2:131-230/310<Extracellular>
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single

coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]

Function:

Odorant receptor.

Subcellular Location:

Cell membrane.

Post-translational modifications:

Belongs to the G-protein coupled receptor 1 family.

SWISS:

Q8NH69

Gene ID:

390148

Database links:

[Entrez Gene: 390148](#) Human

[SwissProt: Q8NH69](#) Human

[Unigene: 554526](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.