




Rabbit Anti-MXRA7 antibody

SL19117R

Product Name:	MXRA7
Chinese Name:	基质重塑相关蛋白7抗体
Alias:	FLJ41492; FLJ41492; FLJ46603; HBV PreS1 transactivated protein 1; Matrix remodelling associated 7; Matrix-remodeling-associated protein 7; MXRA7; MXRA7 HUMAN; PS1TP1; TMAP1; Transmembrane anchor protein 1.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Sheep,
Applications:	WB=1:500-2000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	21kDa
Cellular localization:	The cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human MXRA7:121-204/204
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	

Picture:



Sample: Eye (Mouse) Lysate at 40 ug

Primary: Anti-MXRA7 (SL19117R) at 1/300 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 21 kD

Observed band size: 18 kD