



Rabbit Anti-RHBDD2 antibody

SL19866R

Product Name:	RHBDD2
Chinese Name:	菱形结合域同源蛋白RHBDD2抗体
Alias:	H_RG122E10.2a; H_RG122E10.2b; NPD007; RHBD2_HUMAN; RHBDD2; RHBDDL7; Rhomboid domain containing 2; rhomboid domain containing protein 2; Rhomboid domain-containing protein 2; Rhomboid veinlet like 7; WUGSC:H_RG122E10.2a; WUGSC:H_RG122E10.2b.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Pig,Horse,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	39kDa
Cellular localization:	The cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human RHBDD2:221-320/364
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	