



Rabbit Anti-H5N6 hemagglutinin antibody

SL20019R

Product Name:	H5N6 hemagglutinin
Chinese Name:	A型流感病毒H5N6凝集素
Alias:	hemagglutinin [Influenza A virus (A/mallard/Sweden/40/2002(H5N6))]; Influenza A virus H5N6 hemagglutinin; HA1 (H5N6).
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	H5N6
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	63kDa
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human H5N6 hemagglutinin:351-450/54
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	