



Rabbit Anti-LITD1 antibody

SL20833R

Product Name:	LITD1
Chinese Name:	LITD1蛋白抗体
Alias:	ECAT11; ES cell associated protein 11; ES cell-associated protein 11; FLJ10884; L1td1; LINE 1 type transposase domain containing 1; LINE 1 type transposase domain containing protein 1; LINE-1 type transposase domain-containing protein 1; LITD1_HUMAN; MGC133253; OTTHUMP00000010477; RP5-1155K23.3.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Mouse,Rat,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	99kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human LITD1:651-750/782
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	