



Rabbit Anti-LIX1 antibody

SL2657R

Product Name:	LIX1
Chinese Name:	LIX1抗体
Alias:	C5orf11; Limb expression 1; LIX1; Lix1 homolog (chicken); LIX1_HUMAN; Protein limb expression 1 homolog.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Dog,Cow,Horse,Rabbit,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	32kDa
Cellular localization:	cytoplasmic
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human LIX1:201-282/282
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	The specific function of this protein remains unknown. Similarity: Belongs to the LIX1 family. SWISS:

Q8N485

Gene ID:

167410

Database links:

[Entrez Gene: 167410](#)Human

[Entrez Gene: 66643](#)Mouse

[Omim: 610466](#)Human

[SwissProt: Q8N485](#)Human

[SwissProt: Q6P566](#)Mouse

[Unigene: 656702](#)Human

[Unigene: 268018](#)Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongintech.com