



Rabbit Anti-Neutrophil defensin 4/FITC Conjugated antibody

SL2799R-FITC

Product Name:	Anti-Neutrophil defensin 4/FITC
Chinese Name:	FITC标记的防御素 α 4抗体
Alias:	Neutrophil defensin 4; Corticostatin; Defensin alpha 4; Defensin, alpha 4; HNP-4; HP-4; HNP 4; defensin, alpha 4, corticostatin; neutrophil defensin 4 precursor; neutrophil defensin 4 preproprotein; DEF4; DEFA4; DEF4_HUMAN.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,
Applications:	IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	3.8kDa
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human Defensin alpha 4
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail:	background: Defensins are a family of microbicidal and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. Several alpha defensin genes

are clustered on chromosome 8. This gene differs from other genes of this family by an extra 83-base segment that is apparently the result of a recent duplication within the coding region. The protein encoded by this gene, defensin, alpha 4, is found in the neutrophils; it exhibits corticostatic activity and inhibits corticotropin stimulated corticosterone production. [provided by RefSeq].

Function:

Has antimicrobial activity against Gram-negative bacteria, and to a lesser extent also against Gram-positive bacteria and fungi. Protects blood cells against infection with HIV-1 (in vitro). Inhibits corticotropin (ACTH)-stimulated corticosterone production.

Subunit:

Homodimer.

Subcellular Location:

Secreted.

Similarity:

Belongs to the alpha-defensin family.

Database links:

[Entrez Gene: 1669](#)Human

[Omim: 601157](#)Human

[SwissProt: P12838](#)Human

[Unigene: 591391](#)Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.