



Rabbit Anti-Phospho-ATP1A1 (Ser23) antibody

SL3064R

Product Name:	Phospho-ATP1A1 (Ser23)
Chinese Name:	磷酸化钠钾ATP酶蛋白a1抗体
Alias:	alpha 1 Sodium Potassium ATPase (phospho S23); ATP1A1; Na,K-ATPase alpha-1; Atpa-1; ATPase Na ⁺ /K ⁺ transporting alpha 1 polypeptide; BC010319; EC 3.6.3.9; MGC3285; MGC38419; MGC51750; Na K ATPase alpha A catalytic polypeptide; Na K ATPase catalytic subunit alpha A protein; Na(+)/K(+) ATPase 1; Na(+)/K(+) ATPase alpha-1 subunit; Na ⁺ /K ⁺ ATPase alpha 1 subunit; Na ⁺ /K ⁺ ATPase 1; Nkaa1b; Sodium potassium ATPase alpha 1 polypeptide; Sodium pump 1; Sodium pump subunit alpha-1; Sodium/potassium transporting ATPase alpha 1 chain; Sodium/potassium-transporting ATPase subunit alpha-1.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Mouse,Rat,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	110kDa
Cellular localization:	The cell membrane
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated Synthesised phosphopeptide derived from rat ATP1A1 around the phosphorylation site of Ser23:KK(p-S)KK
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of

	antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	<p>The protein encoded by this gene belongs to the family of P-type cation transport ATPases, and to the subfamily of Na⁺/K⁺-ATPases. Na⁺/K⁺ -ATPase is an integral membrane protein responsible for establishing and maintaining the electrochemical gradients of Na and K ions across the plasma membrane. These gradients are essential for osmoregulation, for sodium-coupled transport of a variety of organic and inorganic molecules, and for electrical excitability of nerve and muscle. This enzyme is composed of two subunits, a large catalytic subunit (alpha) and a smaller glycoprotein subunit (beta). The catalytic subunit of Na⁺/K⁺ -ATPase is encoded by multiple genes. This gene encodes an alpha 1 subunit. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May2009].</p> <p>Function: This is the catalytic component of the active enzyme, which catalyzes the hydrolysis of ATP coupled with the exchange of sodium and potassium ions across the plasma membrane. This action creates the electrochemical gradient of sodium and potassium ions, providing the energy for active transport of various nutrients.</p> <p>Subunit: Interacts with SIK1. Composed of three subunits: alpha (catalytic), beta and gamma. Binds the HLA class II histocompatibility antigen, DR1.</p> <p>Subcellular Location: Cell membrane. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.</p> <p>Tissue Specificity: Belongs to the cation transport ATPase (P-type) (TC 3.A.3) family. Type IIC subfamily.</p> <p>Post-translational modifications: Phosphorylation on Tyr-10 modulates pumping activity. Dephosphorylation by protein phosphatase 2A (PP2A) following increases in intracellular sodium, leading to increase catalytic activity.</p> <p>Similarity: Belongs to the cation transport ATPase (P-type) (TC 3.A.3) family. Type IIC subfamily.</p> <p>SWISS: P06685</p> <p>Gene ID: 24211</p> <p>Database links:</p>

[Entrez Gene: 476](#) Human

[Entrez Gene: 11928](#) Mouse

[Entrez Gene: 24211](#) Rat

[Omim: 182310](#) Human

[SwissProt: P05023](#) Human

[SwissProt: Q8VDN2](#) Mouse

[SwissProt: P06685](#) Rat

[Unigene: 371889](#) Human

[Unigene: 280103](#) Mouse

[Unigene: 217534](#) Rat


[Unigene: 2992](#) Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Channel protein (Channel Protein) 钠钾ATP酶是位于The cell membrane上的一种glycoprotein,与ATP的分解和细胞内外钠、钾离子的转运密切相关,哺乳动物各种组织细胞的钠钾ATP酶的immunology特性基本相同。

Picture:



Sample:

293T(Human) Cell Lysate at 40 ug

Primary: Anti-Phospho-ATP1A1 (Ser23) (SL3064R) at 1/300 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 110 kD

Observed band size: 110 kD