



Rabbit Anti-SMC1 antibody

SL3568R

Product Name:	SMC1
Chinese Name:	染色体结构维持蛋白1抗体
Alias:	CDLS2; SMC1; Segregation of mitotic chromosomes 1; Segregation of mitotic chromosomes like 1; SMC 1; SMC1A; SMC1alpha; SMC1alpha protein; DXS423E; KIAA0178; SB1.8; SMC1B; SMC1BETA; SMC1beta protein; SMC1L1; SMC1L2; SMCB; Structural maintenance of chromosome 1 like 1 protein; Structural maintenance of chromosome 1 like 2 protein; structural maintenance of chromosomes 1-like 2 (yeast); Structural maintenance of chromosomes 1A; Structural maintenance of chromosomes 1B; Structural Maintenance of Chromosomes-1 Like 1; SMC1B_HUMAN.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Chicken,Dog,Pig,Cow,Rabbit,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	136kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human SMC1:151-250/1233
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	Proper cohesion of sister chromatids is a prerequisite for the correct segregation of

chromosomes during cell division. The cohesin multiprotein complex is required for sister chromatid cohesion. This complex is composed partly of two structural maintenance of chromosomes (SMC) proteins, SMC3 and either SMC1B or the protein encoded by this gene. Most of the cohesin complexes dissociate from the chromosomes before mitosis, although those complexes at the kinetochore remain. Therefore, the encoded protein is thought to be an important part of functional kinetochores. In addition, this protein interacts with BRCA1 and is phosphorylated by ATM, indicating a potential role for this protein in DNA repair. This gene, which belongs to the SMC gene family, is located in an area of the X-chromosome that escapes X inactivation. Mutations in this gene result in Cornelia de Lange syndrome. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2013]

Function:

Involved in chromosome cohesion during cell cycle and in DNA repair. Central component of cohesin complex. The cohesion complex is required for the cohesion of sister chromatids after DNA replication. The cohesin complex apparently forms a large proteinaceous ring within which sister chromatids can be trapped. At anaphase, the complex is cleaved and dissociates from chromatin, allowing sister chromatids to segregate. The cohesin complex may also play a role in spindle pole assembly during mitosis. Involved in DNA repair via its interaction with BRCA1 and its related phosphorylation by ATM, or via its phosphorylation by ATR. Works as a downstream effector both in the ATM/NBS1 branch and in the ATR/MSH2 branch of S-phase checkpoint.

Subunit:

Interacts with POLE. Interacts with SYCP2. Interacts with BRCA1. Found in a complex with CDCA5, SMC3 and RAD21, PDS5A/APRIN and PDS5B/SCC-112 (By similarity). Forms a heterodimer with SMC3 in cohesin complexes. Cohesin complexes are composed of the SMC1 (SMC1A or SMC1B) and SMC3 heterodimer attached via their hinge domain, RAD21 which link them, and one STAG protein (STAG1, STAG2 or STAG3), which interacts with RAD21. In germ cell cohesion complexes, SMC1A is mutually exclusive with SMC1B. Interacts with BRCA1. Interacts with NDC80.

Subcellular Location:

Nucleus. Chromosome. Chromosome, centromere, kinetochore. Note=Associates with chromatin. Before prophase it is scattered along chromosome arms. During prophase, most of cohesion complexes dissociate from chromatin probably because of phosphorylation by PLK, except at centromeres, where cohesion complexes remain. At anaphase, the RAD21 subunit of the cohesion complex is cleaved, leading to the dissociation of the complex from chromosomes, allowing chromosome separation. In germ cells, cohesion complex dissociates from chromatin at prophase I, and may be replaced by a meiosis-specific cohesin complex. The phosphorylated form on Ser-957 and Ser-966 associates with chromatin during G1/S/G2 phases but not during M phase, suggesting that phosphorylation does not regulate cohesin function. Integral component of the functional centromere-kinetochore complex at the kinetochore region during mitosis.

Post-translational modifications:

Phosphorylated by ATM upon ionizing radiation in a NBS1-dependent manner.
Phosphorylated by ATR upon DNA methylation in a MSH2/MSH6-dependent manner.
Phosphorylation of Ser-957 and Ser-966 activates it and is required for S-phase checkpoint activation.

DISEASE:

Defects in SMC1A are the cause of Cornelia de Lange syndrome type 2 (CDLS2) [MIM:300590]; also known as Cornelia de Lange syndrome X-linked. CDLS is a clinically heterogeneous developmental disorder associated with malformations affecting multiple systems. CDLS is characterized by facial dysmorphisms, abnormal hands and feet, growth delay, cognitive retardation and various other malformations including gastroesophageal dysfunction and cardiac, ophthalmologic and genitourinary anomalies.

Similarity:

Belongs to the SMC family. SMC1 subfamily.

SWISS:

Q14683

Gene ID:

8243

Database links:

[Entrez Gene: 27127](#)Human

[Entrez Gene: 8243](#)Human

[Entrez Gene: 140557](#)Mouse

[Entrez Gene: 24061](#)Mouse

[Entrez Gene: 300121](#)Rat

[Entrez Gene: 63996](#)Rat

[Oimim: 300040](#)Human

[Oimim: 608685](#)Human

[SwissProt: Q14683](#)Human

[SwissProt: Q8NDV3](#)Human

[SwissProt: Q920F6](#)Mouse

[SwissProt: Q9CU62](#)Mouse

[SwissProt: Q9Z1M9](#)Rat

[Unigene: 334176](#)Human


[Unigene: 182737](#)Mouse

[Unigene: 56972](#)Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

Picture:



Key	Name	Parameter	Gate
—	Hela-blank-20150604.021	FL1-H	G1
—	bs-0295G-FITC-Hela-5.036	FL1-H	G1
—	bs-0295P-(FITC)-Hela-5.037	FL1-H	G1
—	bs-3568R-(FITC)-Hela-5.049	FL1-H	G1

Positive control: HeLa cells

Concentration: 5 μ g/10⁶ cells

Incubation conditions: Avoid light , 30 minutes on the ice.