



## Rabbit Anti-Phospho-eNOS (Thr495) antibody

SL3731R

<b>Product Name:</b>	Phospho-eNOS (Thr495)
<b>Chinese Name:</b>	磷酸化一氧化氮合成酶3(内皮型)抗体
<b>Alias:</b>	eNOS (phospho T494)mo, rat; eNOS (phospho T495)human; p-eNOS (phospho T495)human; eNOS (Phospho-Thr495)human; cNOS; eNOS; Constitutive NOS; EC NOS; ecNOS; Endothelial nitric oxidase synthase; Endothelial nitric oxide synthase; Endothelial nitric oxide synthase 3; Endothelial NOS; Nitric oxide synthase 3 (endothelial cell); Nitric oxide synthase 3; Nitric oxide synthase 3 endothelial cell; Nitric oxide synthase endothelial; nitric oxide synthase, endothelial; NOS 3; NOS III; NOS type III; NOS3; NOSIII; NOS3 HUMAN.
文献引用 PubMed :	<p><b>Specific References(1)</b> SL3731R has been referenced in 1 publications.</p> <p><b>[IF=4.99]</b>Papinska, A. M., et al. "Angiotensin-(1-7) administration benefits cardiac, renal, and progenitor cell function in db/db mice." British Journal of Pharmacology(2015).IHC-P;Mouse.</p> <p style="text-align: right;"><a href="#">PubMed:26075703</a></p>
<b>Organism Species:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>React Species:</b>	Human,Mouse,Rat,Dog,Pig,Cow,Rabbit,Sheep,Guinea Pig,
<b>Applications:</b>	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight:</b>	133kDa
<b>Cellular localization:</b>	cytoplasmicThe cell membrane
<b>Form:</b>	Lyophilized or Liquid
<b>Concentration:</b>	1mg/ml
<b>immunogen:</b>	KLH conjugated Synthesised phosphopeptide derived from human eNOS around the

	phosphorylation site of Thr495:KK(p-T)FK
<b>Lsotype:</b>	IgG
<b>Purification:</b>	affinity purified by Protein A
<b>Storage Buffer:</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Storage:</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>PubMed:</b>	<a href="#">PubMed</a>
<b>Product Detail:</b>	<p>eNOS is a calcium/calmodulin dependent enzyme which undergoes several post translational modifications, including acylation with myristate and palmitate, and phosphorylation on numerous residues. As with the other members of the NOS family, eNOS derives the diffusible multifunctional second messenger NO from L arginine through a series of reactions in which L citrulline is an intermediate. eNOS plays an important role in controlling vascular tone, platelet aggregation, and cardiac myocyte function.</p> <p><b>Function:</b>          Produces nitric oxide (NO) which is implicated in vascular smooth muscle relaxation through a cGMP-mediated signal transduction pathway. NO mediates vascular endothelial growth factor (VEGF)-induced angiogenesis in coronary vessels and promotes blood clotting through the activation of platelets.          Isoform eNOS13C: Lacks eNOS activity, dominant-negative form that may down-regulate eNOS activity by forming heterodimers with isoform 1.</p> <p><b>Subunit:</b>          Homodimer. Interacts with NOSIP and NOSTRIN.</p> <p><b>Subcellular Location:</b>          Cell membrane. Membrane, caveola. Cytoplasm, cytoskeleton. Golgi apparatus.          Note=Specifically associates with actin cytoskeleton in the G2 phase of the cell cycle; which is favored by interaction with NOSIP and results in a reduced enzymatic activity.</p> <p><b>Tissue Specificity:</b>          Platelets, placenta, liver and kidney.</p> <p><b>Post-translational modifications:</b>          Phosphorylation by AMPK at Ser-1177 in the presence of Ca(2+)-calmodulin (CaM) activates activity. In absence of Ca(2+)-calmodulin, AMPK also phosphorylates Thr-495, resulting in inhibition of activity. Phosphorylation of Ser-114 by CDK5 reduces activity.</p> <p><b>Similarity:</b>          Belongs to the NOS family.          Contains 1 FAD-binding FR-type domain.          Contains 1 flavodoxin-like domain.</p>

**SWISS:**  
P29474

**Gene ID:**  
4846

**Database links:**

[Entrez Gene: 287024](#)Cow

[Entrez Gene: 403784](#)Dog

[Entrez Gene: 4846](#)Human

[Entrez Gene: 18127](#)Mouse

[Entrez Gene: 397557](#)Pig

[Entrez Gene: 24600](#)Rat

[Olim: 163729](#)Human


[SwissProt: P29473](#)Cow

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

催化生物体内一氧化氮(NO)生成的酶。分神经型一氧化氮合成的酶(nNOS or NOS-1)、诱导型一氧化氮合成的酶(iNOS or NOS-2)、内皮型一氧化氮合成的酶(eNOS or NOS-3)。此抗体识别分子量为135kDa的内皮型一氧化氮合成的酶。

Picture:



Sample:


K562(Human) Cell Lysate at 30 ug

Primary: Anti-Phospho-eNOS (Thr495) (Ser319) (SL3731R) at 1/300 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 133 kD

Observed band size: 133 kD



Sample:

Lane1: Liver (Mouse) Lysate at 30 ug

Lane2: Kidney (Mouse) Cell Lysate at 30 ug

Primary: Anti-Phospho-eNOS (SL3731R) at 1:300 dilution;

Secondary: HRP conjugated Goat-Anti-Rabbit IgG(bse-0295G) at 1: 3000 dilution;

Predicted band size : 133kD

Observed band size : 135kD