



Rabbit Anti-phospholamban antibody

SL4197R

Product Name:	phospholamban
Chinese Name:	受磷蛋白/心脏磷蛋白抗体
Alias:	Cardiac phospholamban; CMD1P; PLB; PLN; PPLA_HUMAN.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Pig,Cow,
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	5.7kDa
Cellular localization:	cytoplasmicThe cell membrane Mitochondrion
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human phospholamban:1-35/52
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	The Sarco(endo)plasmic-reticulum (SER) regulatory protein, Phospholamban (PLB), is a small, plasma membrane-associated phospho-protein found in the SER of cardiac, smooth and slow-twitch muscle. Believed to assemble into a pentamer, PLB regulates cardiac contractility and Ca ²⁺ affinity for cardiac SER Ca ²⁺ ATPase (SERCA2a). Non-phosphorylated PLB associates with SERCA2a, and inhibits Ca ²⁺ reuptake into the SER. PLB activation occurs when key Serine/Threonine residues in PLB (Ser-10, Ser-16, Thr-17) are phosphorylated by numerous effectors, which include PKC, PKA,

PKG, and CaM kinase. Phosphorylation of PLB causes dissociation from SERCA2a and a subsequent increase in the rate of Ca²⁺ reuptake into the SER, which accelerates ventricular relaxation.

Function:

Phospholamban has been postulated to regulate the activity of the calcium pump of cardiac sarcoplasmic reticulum.

Subcellular Location:

Mitochondrion membrane; Single-pass membrane protein. Sarcoplasmic reticulum.

Tissue Specificity:

Heart.

Post-translational modifications:

Phosphorylated at Thr-17 by CaMK2, and in response to beta-adrenergic stimulation. Phosphorylation by DMPK may stimulate sarcoplasmic reticulum calcium uptake in cardiomyocytes.

DISEASE:

Defects in PLN are the cause of cardiomyopathy dilated type 1P (CMD1P) [MIM:609909]. Dilated cardiomyopathy is a disorder characterized by ventricular dilation and impaired systolic function, resulting in congestive heart failure and arrhythmia. Patients are at risk of premature death.

Defects in PLN are the cause of familial hypertrophic cardiomyopathy type 18 (CMH18) [MIM:613874]. CMH18 is a hereditary heart disorder characterized by ventricular hypertrophy, which is usually asymmetric and often involves the interventricular septum. The symptoms include dyspnea, syncope, collapse, palpitations, and chest pain. They can be readily provoked by exercise. The disorder has inter- and intrafamilial variability ranging from benign to malignant forms with high risk of cardiac failure and sudden cardiac death.

Similarity:

Belongs to the phospholamban family.

SWISS:

P26678

Gene ID:

5350

Database links:

[Entrez Gene: 5350](#)Human

[Entrez Gene: 18821](#)Mouse

[Entrez Gene: 64672](#)Rat

[Omim: 172405](#)Human

[SwissProt: A4IFH6](#)Cow

[SwissProt: P26678](#)Human

[SwissProt: P61014](#)Mouse

[SwissProt: P61013](#)Pig

[SwissProt: P61015](#)Rabbit

[SwissProt: P61016](#)Rat

[Unigene: 170839](#)Human

[Unigene: 34145](#)Mouse

[Unigene: 9740](#)Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com