



Rabbit Anti-ZNF641 antibody

SL4342R

Product Name:	ZNF641
Chinese Name:	Zinc finger protein641抗体
Alias:	DKFZp667D1012; FLJ31295; RGD1310412; Zfp641; Zinc finger protein 641.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Dog,Pig,Cow,Horse,Rabbit,Sheep,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	50kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human ZNF641:81-180/438
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	ZNF641 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 5 C2H2-type zinc fingers and 1 KRAB domain. ZNF641 activates transcriptional activities of SRE and AP-1. Subcellular Location: Nucleus.

SWISS:
Q96N77

Gene ID:
121274

Database links:

[Entrez Gene: 538714](#) Cow

[Entrez Gene: 121274](#) Human

[Entrez Gene: 239652](#) Mouse

[Entrez Gene: 300197](#) Rat

[SwissProt: Q96N77](#) Human


[SwissProt: Q8BZ34](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com

Picture:



Sample:

HL-60 (Human) Cell Lysate at 30 ug

Primary: Anti-ZNF641 (SL4342R) at 1/1000 dilution

Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution

Predicted band size: 50 kD

Observed band size: 50 kD

www.sunlongbiotech.com