



Rabbit Anti-ZNF660 antibody

SL4350R

Product Name:	ZNF660
Chinese Name:	Zinc finger protein660抗体
Alias:	Zinc finger protein 660.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	38kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human ZNF660:1-100/331
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	<p>ZNF660 belongs to the krueppel C2H2-type zinc-finger protein family. It contains ten C2H2-type zinc fingers. ZNF660 may be involved in transcriptional regulation.</p> <p>Subcellular Location: Nucleus.</p> <p>SWISS:</p>

Q6AZW8

Gene ID:
285349

Database links:

[Entrez Gene: 285349](#) Human

[SwissProt: Q6AZW8](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com