



Rabbit Anti-ZNF710 antibody

SL4373R

Product Name:	ZNF710
Chinese Name:	Zinc finger protein710抗体
Alias:	DKFZp547K1113; FLJ00306; FLJ37393; MGC163413; Zinc finger protein 710; ZN710 HUMAN; Znf710.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Dog,Pig,Cow,Horse,Sheep,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	75kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human ZNF710:401-500/664
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	ZNF710 (Zinc Finger Protein 710) is a Protein Coding gene. Among its related pathways are Gene Expression. GO annotations related to this gene include nucleic acid binding and transcription cofactor activity. An important paralog of this gene is ZNF366. Function:

May be involved in transcriptional regulation.

Subcellular Location:

Nucleus.

Similarity:

Belongs to the krueppel C2H2-type zinc-finger protein family.
Contains 11 C2H2-type zinc fingers.

SWISS:

8N1W2

Gene ID:

374655

Database links:

[Entrez Gene: 374655](#) Human

[Entrez Gene: 209225](#) Mouse

[Entrez Gene: 293044](#) Rat

[GenBank: NM_198526](#) Human

[SwissProt: Q8N1W2](#) Human

[SwissProt: Q3U288](#) Mouse

[Unigene: 459311](#) Human

[Unigene: 290925](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.