



Rabbit Anti-PVRIG antibody

SL4439R

| | |
|-------------------------------|---|
| Product Name: | PVRIG |
| Chinese Name: | 脊髓灰质炎病毒受体相关结构域蛋白抗体 |
| Alias: | Poliovirus receptor-related immunoglobulin domain-containing protein; PVRIG; C7orf15; PVRIG_HUMAN; Transmembrane protein PVRIG. |
| Organism Species: | Rabbit |
| Clonality: | Polyclonal |
| React Species: | Human, |
| Applications: | ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Molecular weight: | 34kDa |
| Cellular localization: | The cell membrane |
| Form: | Lyophilized or Liquid |
| Concentration: | 1mg/ml |
| immunogen: | KLH conjugated synthetic peptide derived from human PVRIG:151-250/326 |
| Lsotype: | IgG |
| Purification: | affinity purified by Protein A |
| Storage Buffer: | Preservative: 15mM Sodium Azide, Constituents: 1% BSA, 0.01M PBS, pH 7.4 |
| Storage: | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| PubMed: | PubMed |
| Product Detail: | |