



## Rabbit Anti-LOC103437570 antibody

SL4517R

|                          |   |
|--------------------------|---|
| <b>Product Name:</b>     | LOC103437570  |
| <b>Chinese Name:</b>     | LOC103437570抗体  |
| <b>Alias:</b>            | ATP-dependent Clp protease ATP-binding subunit ClpA homolog CD4B, chloroplastic; ATP-dependent Clp protease ATP-binding subunit ClpA homolog CD4B, chloroplastic [ Malus domestica (apple) ]; Apple ClpA/B.   |
| <b>Organism Species:</b> | Rabbit  |
| <b>Clonality:</b>        | Polyclonal  |
| <b>React Species:</b>    | Arabidopsis Thaliana,Apple  |
| <b>Applications:</b>     | ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800 (Paraffin sections need antigen repair)<br>not yet tested in other applications.<br>optimal dilutions/concentrations should be determined by the end user.   |
| <b>Molecular weight:</b> | 103kDa  |
| <b>Form:</b>             | Lyophilized or Liquid   |
| <b>Concentration:</b>    | 1mg/ml  |
| <b>immunogen:</b>        | KLH conjugated synthetic peptide derived from Apple LOC103437570:   |
| <b>Lsotype:</b>          | IgG   |
| <b>Purification:</b>     | affinity purified by Protein A  |
| <b>Storage Buffer:</b>   | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.  |
| <b>Storage:</b>          | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| <b>PubMed:</b>           | <a href="#">PubMed</a>  |
| <b>Product Detail:</b>   |   |