



Rabbit Anti-Duck salmonella antibody

SL4805R

Product Name:	Duck salmonella
Chinese Name:	鸭沙门氏菌抗体
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	salnella
Applications:	ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	Whole cell protein of Duck salmonella:
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	鸭沙门氏菌病又称鸭副伤寒, 是由沙门氏菌属的细菌引起鸭的急性或慢性传染病。其特点是污染面大, 雏鸭特别易感, 多呈急性经过, 而成年鸭多呈慢性或隐性感染成为带菌者。