




Rabbit Anti-CHCHD2 antibody

SL6407R

Product Name:	CHCHD2
Chinese Name:	丙型肝炎NS2反式调节蛋白/衰老相关的基因10蛋白抗体
Alias:	mitochondrial; 16.7kD protein; Aging Associated Gene 10 Protein; Aging-associated gene 10 protein; C7orf17; CHCH2_HUMAN; CHCHD 2; CHCHD2; Coiled Coil Helix Coiled Coil Helix Domain Containing 2; Coiled coil helix coiled coil helix domain containing protein 2, mitochondrial; Coiled-coil-helix-coiled-coil-helix domain-containing protein 2; HCV NS2 trans regulated protein; HCV NS2 trans-regulated protein; NS2TP.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,Chicken,Dog,Pig,Horse,Rabbit,Guinea Pig,G
Applications:	IHC-P=1:100-200IHC-F=1:100-200 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	12 kDa
Cellular localization:	cytoplasmic Mitochondrion
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human CHCHD2:96-130/151
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	



Picture:

Paraformaldehyde-fixed, paraffin embedded (Mouse stomach); Antigen retrieval by boiling in sodium citrate buffer (pH6.0) for 15min; Block endogenous peroxidase by 3% hydrogen peroxide for 20 minutes; Blocking buffer (normal goat serum) at 37°C for 30min; Antibody incubation with (CHCHD2) Polyclonal Antibody, Unconjugated (SL6407R) at 1:400 overnight at 4°C, followed by a conjugated secondary antibody (sp-0023) for 20 minutes and DAB staining.