



## Rabbit Anti-ZNF121 antibody

SL7148R

<b>Product Name:</b>	ZNF121
<b>Chinese Name:</b>	Zinc finger protein121抗体
<b>Alias:</b>	D19S204; ZHC 32; ZHC32; Zinc finger protein 121; Zinc finger protein 20; ZNF 121; ZNF 20; ZNF20.
<b>Organism Species:</b>	Rabbit
<b>Clonality:</b>	Polyclonal
<b>React Species:</b>	Human,
<b>Applications:</b>	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight:</b>	45kDa
<b>Cellular localization:</b>	The nucleus
<b>Form:</b>	Lyophilized or Liquid
<b>Concentration:</b>	1mg/ml
<b>immunogen:</b>	KLH conjugated synthetic peptide derived from human ZNF121:21-120/390
<b>Lsotype:</b>	IgG
<b>Purification:</b>	affinity purified by Protein A
<b>Storage Buffer:</b>	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
<b>Storage:</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>PubMed:</b>	<a href="#">PubMed</a>
<b>Product Detail:</b>	ZNF121 belongs to the krueppel C2H2-type zinc-finger protein family. It contains 11 C2H2-type zinc fingers. ZNF121 may be involved in transcriptional regulation. <b>Function:</b> May be involved in transcriptional regulation.

**Subcellular Location:**

Nuclear.

**Similarity:**

Belongs to the krueppel C2H2-type zinc-finger protein family.  
Contains 11 C2H2-type zinc fingers.

**SWISS:**

P58317

**Gene ID:**

7675

**Database links:**

[Entrez Gene: 7675](#) Human

[Omim: 194628](#) Human

[SwissProt: P58317](#) Human

[Unigene: 501537](#) Human

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlongbiotech.com