



Rabbit Anti-ZNF251 antibody

SL7196R

Product Name:	ZNF251
Chinese Name:	Zinc finger protein251抗体
Alias:	Zinc finger protein 251; Zfp251; ZNF251.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800ICC=1:100-500IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	71kDa
Cellular localization:	The nucleus
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human ZNF251:101-200/671
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	May be involved in transcriptional regulation. Subcellular Location: Nuclear Similarity: Belongs to the krueppel C2H2-type zinc-finger protein family.

Contains 14 C2H2-type zinc fingers.
Contains 1 KRAB domain.

SWISS:
Q9BRH9

Gene ID:
90987

Database links:

[Entrez Gene: 90987](#) Human

[Entrez Gene: 71591](#) Mouse

[Entrez Gene: 366954](#) Rat

[SwissProt: Q9BRH9](#) Human

[SwissProt: Q6PCX8](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

www.sunlabiotech.com