



Rabbit Anti-KLHL26 antibody

SL8052R

Product Name:	KLHL26
Chinese Name:	Kelch样蛋白26抗体
Alias:	FLJ11078; Kelch like 26 (Drosophila); Kelch like protein 26.
Organism Species:	Rabbit
Clonality:	Polyclonal
React Species:	Human,Mouse,Rat,
Applications:	WB=1:500-2000ELISA=1:500-1000IHC-P=1:400-800IHC-F=1:400-800IF=1:100-500 (Paraffin sections need antigen repair) not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight:	68kDa
Cellular localization:	The nucleuscytoplasmicThe cell membraneExtracellular matrixSecretory protein
Form:	Lyophilized or Liquid
Concentration:	1mg/ml
immunogen:	KLH conjugated synthetic peptide derived from human KLHL26:131-240/615
Lsotype:	IgG
Purification:	affinity purified by Protein A
Storage Buffer:	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage:	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
PubMed:	PubMed
Product Detail:	KLHL26 has three named isoforms. The function of KLHL26 remains unknown. SWISS: 55295 Gene ID: Q53HC5

Database links:

[Entrez Gene: 55295](#)Human

[Entrez Gene: 234378](#)Mouse

[Entrez Gene: 290657](#)Rat

[SwissProt: Q53HC5](#)Human

[SwissProt: Q8BGY4](#)Mouse

[SwissProt: Q5RGB8](#)Zebrafish

[Unigene: 250632](#)Human

[Unigene: 187090](#)Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.