



## Rabbit Anti-CASD1/C7orf12 antibody

SL9403R

|                               |   |
|-------------------------------|---|
| <b>Product Name:</b>          | CASD1/C7orf12   |
| <b>Chinese Name:</b>          | 第七号染色体开放阅读框12抗体   |
| <b>Alias:</b>                 | C7orf12; CAS1 domain-containing protein 1; Casd1; CASD1_HUMAN; FLJ21213; FLJ21879; FLJ41901; NBLA04196.   |
| <b>Organism Species:</b>      | Rabbit  |
| <b>Clonality:</b>             | Polyclonal  |
| <b>React Species:</b>         | Human,Mouse,Rat,Chicken,Dog,Pig,Cow,Horse,Sheep,  |
| <b>Applications:</b>          | ELISA=1:500-1000IHC-F=1:400-800IF=1:50-200 (Paraffin sections need antigen repair)<br>not yet tested in other applications.<br>optimal dilutions/concentrations should be determined by the end user.   |
| <b>Molecular weight:</b>      | 88kDa   |
| <b>Cellular localization:</b> | The cell membrane   |
| <b>Form:</b>                  | Lyophilized or Liquid   |
| <b>Concentration:</b>         | 1mg/ml  |
| <b>immunogen:</b>             | KLH conjugated synthetic peptide derived from human CASD1/C7orf12:201-300/797   |
| <b>Lsotype:</b>               | IgG   |
| <b>Purification:</b>          | affinity purified by Protein A  |
| <b>Storage Buffer:</b>        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.  |
| <b>Storage:</b>               | Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| <b>PubMed:</b>                | <a href="#">PubMed</a>  |
| <b>Product Detail:</b>        |   |