



SunLong Biotech Co.,LTD
Tel: 0086-571- 56623320 Fax:0086-571- 56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Mouse Anti-CD73 antibody

Catalog Number	SLM-51239M
Size	100ul
Concentration	1mg/ml
Host	Mouse
Target Protein	CD73
Clonality	Monoclonal
Isotype	Mouse IgG1
Entrez Gene	4907
Swiss Prot	P21589
Source	Recombinant human CD73 protein
Purification	affinity purified by Protein A
Storage	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles..
Background	The protein encoded by this gene is a plasma membrane protein that catalyzes the conversion of extracellular nucleotides to membrane-permeable nucleosides. The encoded protein is used as a determinant of lymphocyte differentiation. Defects in this gene can lead to the calcification of joints and arteries. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Mar 2011]
Applications	WB=1:300-1000,IHC-P=1:200-400,Flow-Cyt=1ug/Test
Cross Reactive Species	Human

For research use only. Not intended for diagnostic or therapeutic use.